

# Virtually Particle-Free Rt<sup>®</sup>-Silica BOND Columns Provide Reliable PLOT Column Performance With Less Time Lost for Maintenance

By Amanda Rigdon, Bill Bromps, Tom Veza, and Jaap de Zeeuw, Restek  
Restek Corporation, 110 Benner Circle, Bellefonte, PA 16823 USA  
Email: amanda.rigdon@restek.com • Web: www.restek.com • Tel: 1-814-353-1300 ext. 2422

- Optimized manufacturing process practically eliminates particle release, reducing downtime due to system obstructions and damage from particles.
- Bonded silica stationary phase minimizes impact of water, resulting in reproducible retention times for water-containing samples.
- Versatile, highly retentive column ideal for analysis of light hydrocarbons, sulfur gases, halocarbons, and carbon dioxide at temperatures above ambient.
- Individually QC tested with sensitive unsaturated C4 probes to ensure consistent selectivity

Porous layer open tubular (PLOT) columns are very useful to GC analysts working on a wide variety of applications, and the unique selectivity of PLOT columns makes them particularly good for separating gaseous compounds without cryogenic cooling. However, the overall utility of traditional PLOT columns is hampered by the characteristic instability of the porous layer that coats the inside of the column. With most PLOT columns, particles that shed from the porous layer create significant problems because they can form obstructions inside the column that can alter flow, causing retention time instability. In addition, particle build-up makes frequent maintenance necessary as jets become obstructed and detectors become contaminated. In contrast, new Rt<sup>®</sup>-Silica BOND columns from Restek are exceptionally robust due to optimized manufacturing and phase bonding steps that practically eliminate particle release. This exceptional stability—in combination with high loadability, inertness, and consistent selectivity—makes these new columns extremely reliable and ideal for the analysis of light hydrocarbons, sulfur gases, and halocarbons. In addition, carbon dioxide and other permanent gases can be retained at ambient temperature on this silica-based column. This article demonstrates the robustness of the Rt<sup>®</sup>-Silica BOND column and its performance for many of the applications relevant to testing natural gas and light hydrocarbon streams.

## Virtually Particle-Free and Water Resistant PLOT Performance

Restek's proprietary manufacturing technique for the Rt<sup>®</sup>-Silica BOND column results in an extremely stable porous layer with traditional PLOT column loadability and retention without loose particles that can damage valves and foul FID jets.

Figure 1 shows a magnified picture of three fused silica columns. The middle column is a traditionally-manufactured PLOT column, the bottom column is a wall-coated open tubular column, and the top column is an Rt<sup>®</sup>-Silica BOND PLOT column. Note the uneven layer of particles on the middle column, as well as areas where the particles have completely detached from the column wall; this causes irregularities in the internal diameter of the column that can cause retention time instability. In comparison, the Rt<sup>®</sup>-Silica BOND column looks identical to the wall-coated open tubular column, with no visible shedding of particles or peeling of the coating layer. While the Rt<sup>®</sup>-Silica BOND column does contain a porous layer, the structure of this layer is extremely fine and well-adhered to the column wall, ensuring virtually particle-free operation over the lifetime of the column.





Figure 1: Traditional PLOT columns (middle) have an uneven coating of particles that can shed, fouling instrument parts. Rt<sup>®</sup>-Silica BOND columns (top) have a very fine porous layer with no visible particles and look very similar to wall-coated open tubular columns (bottom).

The manufacturing process used to make Rt<sup>®</sup>-Silica BOND columns results in a PLOT column with high selectivity, retention, and capacity without the particle shedding associated with conventional PLOT columns. This provides improved column robustness and less downtime for maintenance. The particle-free nature of this column is evidenced by a particle-generation experiment in which a column was temperature and pressure ramped multiple times. Changes in temperature cause changes in pressure, which result in particle shedding in traditional PLOT columns. Free particles generate large spikes when they hit the flame ionization detector (FID), interfering with quantification. In addition, the particles themselves can obstruct FID jets and damage valves. Note that no large particle spikes were generated when this experiment was carried out on a brand new Rt<sup>®</sup>-Silica BOND column (Figure 2).



<b>Column</b>	Rt <sup>®</sup> -Silica BOND, 30 m, 0.32 mm ID (cat.# 19785)
<b>Injection</b>	split (split ratio 35:1)
<b>Liner:</b>	Sky® 2.0 mm ID straight inlet liner (cat.# 23313.1)
<b>Inj. Temp.:</b>	250 °C
<b>Oven</b>	
<b>Oven Temp.:</b>	50 °C to 250 °C at 35 °C/min (hold 5 min) to 50 °C at 70 °C/min
<b>Carrier Gas</b>	He, constant flow
<b>Linear Velocity:</b>	114 cm/sec
<b>Detector</b>	FID @ 260 °C
<b>Make-up Gas</b>	
<b>Flow Rate:</b>	50 mL/min
<b>Make-up Gas</b>	
<b>Type:</b>	N <sub>2</sub>
<b>Hydrogen flow:</b>	4.0 mL/min
<b>Air flow:</b>	4.00 mL/min
<b>Data Rate:</b>	10 Hz
<b>Instrument</b>	Agilent 7890A GC

Figure 2: The Rt<sup>®</sup>-Silica BOND PLOT column shows no large particle spikes, even with temperature and pressure variation.

Another benefit of Restek's proprietary manufacturing process for the Rt®-Silica BOND column is that the stationary phase of the column is composed almost entirely of silica. While silica retains water, it does not adsorb it. Some PLOT materials adsorb water, which changes the retention and selectivity of the column. After analyzing samples containing water, these PLOT columns require extensive thermal conditioning (bakeout) to return their original retention and selectivity. Figure 3 shows a mixture of saturated and unsaturated hydrocarbons analyzed on the Rt®-Silica BOND column both before exposure to water and then immediately after 10 large volume water injections. Even under these experimental conditions of extreme overwetting, the retention and selectivity of the column remain very similar and under normal use conditions would be effectively identical. This consistent water-resistant performance allows analysts to save time by minimizing maintenance and eliminating the extensive bakeout periods associated with other PLOT columns.




Figure 3: Repeated water injections have minimal impact on Rt®-Silica BOND column selectivity and retention, meaning, water-containing samples can be analyzed without requiring time-consuming thermal reconditioning.

## Versatile Column for Many Applications

The new Rt®-Silica BOND column combines the retention, capacity, and selectivity of traditional PLOT columns with virtually particle-free, water-resistant performance. The bonded silica surface provides excellent retention for light hydrocarbons (Figure 4), permanent gases, and halocarbons, allowing for easy analysis of impurities in light hydrocarbon streams. In addition to light hydrocarbon analysis, the Rt®-Silica BOND column is especially selective for sulfur compounds in hydrocarbon streams. Figures 5 and 6 illustrate good separation of sulfur compounds in propane and butane, respectively.




Figure 4: Saturated and unsaturated hydrocarbons are resolved and retained well on the Rt®-Silica BOND column.




Figure 5: Sulfur Compounds in Propane





Figure 6: Sulfur Compounds in Butane

## Rigorous QC Testing Ensures Ultimate Column-to-Column Reproducibility

While column-to-column reproducibility is a must for all analysts, it is especially important in valve, backflushing, or columnswitching applications. With this in mind, a special QC test was designed for the Rt®-Silica BOND column. Performance parameters, including efficiency, selectivity (RI), retention (k), and inertness are evaluated for each and every column. While the QC tests from some manufacturers include some of these parameters, the compounds used to measure RI are not well-retained and not sensitive to changes in column selectivity. The RI compounds used in the QC test for the Rt®-Silica BOND are 1,3-butadiene and methyl acetylene, which are not only very sensitive probes for selectivity, but are of high interest to many analysts. Additionally, while some commercially available PLOT columns are not evaluated for inertness, Rt®-Silica BOND column inertness is measured with propylene, which is a more active, unsaturated hydrocarbon. This QC testing ensures the highest level of column-to-column reproducibility available in the industry for PLOT columns. Figure 7 shows QC results from three separate lots of Rt®-Silica BOND columns.





<b>Peaks</b>	<b>Column</b>	Rt®-Silica BOND, 30 m, 0.32 mm ID (cat.# 19753)	<b>Inj. Temp:</b> 250 °C
1. Methane	<b>Sample</b>	Custom gas standard	<b>Oven Temp:</b> 90 °C (hold 20 min)
2. Propylene	<b>Diluent:</b> Nitrogen		<b>Carrier Gas:</b> H <sub>2</sub> , constant flow
3. n-Butane	<b>Conc.:</b> 1 mole percent each component		<b>Linear Velocity:</b> 28 cm/sec
4. 1,2-Dichlorotetrafluoroethane (CFC-114)	<b>Injection</b>	15 µl, split (split ratio 35:1)	<b>Detector</b>
5. Methylacetylene	<b>Inj. Vol.:</b> Sky® 2.0 mm ID straight inlet liner (cat.# 23313.3)		FID @ 260 °C
6. 1,3-Butadiene			<b>Instrument</b>
7. n-Pentane			Agilent HP6890 GC

Figure 7: Rigorous QC testing ensures column-to-column reproducibility.

The Rt®-Silica BOND column gives you the retention and capacity you need from PLOT columns, along with virtually particle-free and water-resistant operation. The combination of rugged manufacturing and rigorous QC testing ensures every Rt®-Silica BOND column will provide optimal performance and reliable results for every analysis, while minimizing downtime due to maintenance from particle shedding or time-consuming bakeouts due to water contamination. The column's unique selectivity makes it ideal for analysis of hydrocarbons, halogenated compounds, and sulfur gases.

Do you have an article for Petro Industry News? If so, email our editor, Rachael Simpson, today.

[rachael@envirotechpubs.com](mailto:rachael@envirotechpubs.com)

Read, Print, Share or Comment on this Article at: [Petro-Online.com/Articles](http://Petro-Online.com/Articles)



## Enthusiastic Welcome for New Hydrocarbon Testing Show

**PEFTEC** 18th - 19th NOVEMBER 2015 ANTWERP BELGIUM  
 CONFERENCE EXHIBITION & WORKSHOPS  
 PETROLEUM, REFINING & ENVIRONMENTAL MONITORING TECHNOLOGIES

PEFTEC 2015, the international petroleum, refining and environmental event (Antwerp 18th-19th November), is set to become a major date in the diaries of anyone involved in the testing, analysis and monitoring of downstream oil products and petrochemicals. PEFTEC organiser Marcus Pattison says: "The response to our Call for

Papers has been overwhelming and the exhibition is rapidly filling with this sector's major instrumentation and laboratory equipment manufacturers, so we are looking forward to an extremely busy show."

PEFTEC 2015 is a new Conference and Exhibition created specifically for chemists, scientists, laboratory staff, academics, researchers, process engineers and environmental managers working in petrochemicals and in refineries. Taking place at the Antwerp Expo, PEFTEC will be situated at the centre of one the world's largest cluster of petrochemical companies.

Anyone interested in attending should bookmark [www.PEFTEC.com](http://www.PEFTEC.com) for regular updates on the Conference, Seminars and Exhibition.

email: [35067pr@reply-direct.com](mailto:35067pr@reply-direct.com)

## A New Dimension of Spectral Resolution in ICP-OES



Taking pictures with your camera you may find that images often look okay on a first glimpse, but when zooming in fine details are barely noticeable. In contrast, the high information density of high resolution images will reveal even the smallest details.

Since the introduction of ICP optical emission spectrometry, the spectral resolution has been a main concern with respect to accuracy and sensitivity. The wealth of emission lines from the plasma frequently leads to overlapping between lines of the analyte and that of accompanying elements. In many sample matrices such spectral interferences impair the recognition of the analyte signal and thus restrict the free choice of analytical lines. Hence, less sensitive alternative lines are frequently used that exhibit poorer detection limits.

With the introduction of the High-Resolution Array ICP-OES PlasmaQuant PQ 9000, Analytik Jena (Germany) is taking optical emission spectrometry to a new level. Thanks to its High-Resolution Optics, PlasmaQuant PQ 9000 exceeds the spectral resolution of conventional instruments by up to factor 4, which translates into uncompromised accuracy and advanced analytical capabilities as previously obscured analyte lines become detectable.

Along with an unprecedented Automatic Baseline Correction (ABC) protocol and an elaborate software tool for the Correction of Spectral Interferences (CSI) the High-Resolution Optics of PlasmaQuant PQ 9000 brings unconditional confidence to your most delicate analytical routine.

email: [33437pr@reply-direct.com](mailto:33437pr@reply-direct.com)

**L-K Industries** Oil Testing Equipment Supplies Worldwide Since 1930

### Criterion Water Bath System

- Certified for operations in Class 1, Group D, Division 2 areas
- Digitally controlled Temperature (Either degrees F/C)
- Maintains the programmed temperature within specified ASTM/API Temperature Guidelines

6952 Lawndale Street, Houston, TX 77023-25995  
 Tel: (713) 926-2623 • Fax: (713) 926-7736  
 Email: [lkind@lk-ind.com](mailto:lkind@lk-ind.com)

**Web: [www.lk-ind.com](http://www.lk-ind.com)**

email: [4036ad@reply-direct.com](mailto:4036ad@reply-direct.com)