

## Not Your Father's Raman Spectrometer New Handheld Brings Lab Power to Field Petrochemical Applications

James Pasmore, Senior Advisor to the President, Rigaku Corporation

Handheld Raman spectrometers have gained wide acceptance in recent years for identifying unknowns in field applications. They have proven to be especially effective in raw materials inspection for petrochemical, chemical, pharmaceutical, and electronic industries. Raman spectroscopy is a non-invasive, non-destructive analytical technique. Unlike the more popular and closely related technique FTIR, it requires no sample preparation and is insensitive to interference from water or moisture.

In the petroleum industry, these spectrometers help to identify and quantify pure materials and mixtures in liquid and solid forms.

This non-destructive technique offers analytical information previously not available in a handheld system.

The handheld Raman has application to refined oil products, such as methanol content in gasoline, and with a new 1064 nm laser excitation source it is particularly useful for crude and biofuels testing. Engine oil wear analysis has also proven to be an application that the 1064 nm wavelength is most suited to analyse effectively.

As impressive as these applications are, they are qualitative in nature. An example of the use of handheld Raman spectrometer for quantitative work, namely measuring methanol component in blended gasoline, is presented below.

Quantitative results are presented in Figure 4, which covers the entire blending range: from neat gasoline to pure methanol. By choosing appropriate spectral features, a simple linear correlation can be found between methanol percentage and the spectral quantity. With the established equation, the percentage of methanol of known samples is predicted to be within  $\pm 2\%$  of their actual values; for example, 8.3% was predicted for a 10% sample and 51% was predicted for a 50% sample from our own experiments.

In conclusion, although it is widely accepted that handheld Raman can identify a wide range of molecular species, we have demonstrated that an ultra-compact, field deployable Raman spectrometer can also be reliably used to measure quantitatively the methanol content in pre-blended gasoline samples.

TM Xantus and FirstGuard are trademarks of BaySpec Corp

Xantus™



Figure 1. The "experimental setup." The experiment consists of holding the laser power button for a few seconds to shine laser light on the sample vial. In this case is the experiment takes approximately four seconds.

FirstGuard™



Figure 2. An alternate configuration with improved ergonomics for point and shoot applications.

Raman Spectra




Figure 3. Typical Raman spectra taken with Rigaku's Handheld. Integration time: 4 seconds.




Figure 4. Correlation of a spectral quantity with the methanol percentage in blended gasoline. Data were fit with linear regression. Error bars represent  $\pm 3$ .

Reader Reply Card No 83

## Improve control of your process parameters with quick turnaround real time analysis

Drilling rigs, refineries, pipelines and gas pumps

- XRF for elemental analysis of metals and sulfur from ppb to %
- XRT for on-line sulfur analysis from low to high % levels
- Raman for non-destructive qualitative materials identification and quantifiable process elucidation

Robust: Rigorous: Rigaku

Supermini WDXRF multi-element




Micro-Z sulfur™ WDXRF sulfur analyzer



NEX CG EDXRF



FirstGuard™ Raman handheld



NEX XT XRT sulfur gauge



ZSX Petro 3 kW or 4 kW WDXRF



Petro-Pak™: Applications package for petro applications  
Oil-Pak™: Applications for lube oil additives and wear metals

Rigaku

Rigaku Americas Corporation phone: +1 (281) 362-2300 e-mail: info@rigaku.com  
Applied Rigaku Technologies phone: +1 (512) 225-1796 e-mail: info@rigakuedxrf.com  
Rigaku Corporation phone: +81 3-3479-0618 e-mail: rintyo@rigaku.co.jp  
Rigaku Europe phone: +[44] 1732 763367 email: info@rigaku.com

Reader Reply Card No 84

www.rigaku.com/contact/pin.html