



ANALYSIS OF REFINERY GAS STREAMS USING THE SCION INSTRUMENTS REFINERY GAS ANALYSER

Introduction

SCION Instruments offers a solution for refinery gas streams: The Refinery Gas Analyser, RGA. The RGA determines the chemical composition of reformed gas and similar gaseous mixtures containing the following components: Helium, Hydrogen, Oxygen/Argon (not separated), Nitrogen, Carbon monoxide, Carbon dioxide, C1-C6+ Hydrocarbons isomers and Hydrogen Sulphide, in a single analysis under 30 minutes. CompassCDS software will provide complete analyser control, data acquisition and flexible report generation. Moreover, The Rapid RGA provides a substantial reduction in overall analysis time; less than six minutes including H₂S. This is another analyser in the SCION Instruments portfolio.

The RGA configuration is optimised for compliance with ASTM D1946 and ASTM D2163. Other methods which can be met by the SCION Instruments RGA are for example ASTM D1945, ASTM D2504, ASTM D7833, ASTM D2593, ASTM D4424, IP405/ISO7941/EN27941, EN15985 and UOP539, although minor modifications may be required. Additionally, for enhanced sample introduction, the RGA can be optionally equipped with an integrated micro-gasifier. This sample conditioning device ensures complete vaporisations of LPG's and high-pressure samples to prevent any sample discrimination prior to injection. Another option for the RGA is the use of liquid sampling valves. The chemical composition of previous mentioned components can be used to calculate important physical properties of the gas, such as heating (calorific) value, and relative density. Physical properties can be calculated with specialist Gas Calculation Software, Eclipse.

Experimental

The standard RGA is a three channel 456-GC with a multi valve design using both capillary and packed columns. Figure 1 shows the Refinery Gas Analyser. The first channel is optimised for the analysis of helium and hydrogen, with the second channel optimised for the analysis of permanent gases. Finally, the third channel is designed for light hydrocarbons.

For the Helium-Hydrogen channel, the sample is injected by means of a gas-sampling valve onto a Hayesep column. The first fraction containing Helium and Hydrogen is flushed onto a Molecular Sieve column; the rest is backflushed to the vent. Helium and Hydrogen are detected by a TCD.

The Permanent Gas channel is developed for the simultaneous determination of CO₂, the C2 isomers, H₂S, Oxygen/Argon, Nitrogen, Methane and CO. The sample is injected by means of a gas-sampling valve onto a series of Hayesep columns. The fraction containing Oxygen, Nitrogen, CO and Methane is flushed onto a Molecular Sieve column and parked. CO₂, the C2 isomers and H₂S elute to the TCD, bypassing the Molecular Sieve column. After the elution of H₂S, the Molecular Sieve column is set in flow again,

giving the separation of Oxygen/Argon, Nitrogen, Methane and CO.

The Light Hydrocarbon channel determines the range of low boiling hydrocarbons and isomers, C1 – C5 and C6+ (as composite peak). The sample is injected by means of a gas-sampling valve onto the pre-column (SCION-1). The components with a boiling point of below C6 are transferred onto the Al₂O₃/Na₂SO₄ column. The highly selective column separates all individual isomers from the light hydrocarbon fraction. The fraction with a boiling point of C6 and higher is backflushed from the pre-column to the FID, resulting in a peak representing the higher boiling fraction.

The backflush time determines the range of hydrocarbons which will be separated on the Al₂O₃/Na₂SO₄ column.

When liquid sample valves are added to the configuration, LPG's can also be injected onto the RGA.

Results

Figures 2a, 2b and 2c show the chromatograms obtained from the three RGA channels, when a Refinery Gas Standard was analysed.

CompassCDS software provides complete analyser control, data acquisition and flexible reporting. The composition of the sample will give important information about the physical properties. To calculate the physical properties, Eclipse Software was used.

CompassCDS was optimised to integrate Eclipse results automatically, generating a physical properties report. The physical properties of interest can be selected from a list for full customisation and reporting.

The software uses specific peak identification tables to determine the components on each of the three channels. Every channel has its own peak identification table. Each peak in this table can be calibrated on detector response using a single or multi-level calibration technique.

The peaks of a sample are identified, and the concentrations are calculated using the calibration data. All peaks can be normalised to 100 mol%. The normalised results are used to calculate the




Figure 1: SCION 456 Refinery Gas Analyser.

desired physical properties.


For the analysis of refinery gas, application EN15984 is used as an example. Figure 3 profiles a typical report generated, customise to report all physical properties possible for that application.

Eclipse has also a special enhancement for oxygen/nitrogen correction. It recalculates the concentration of nitrogen based on the amount of detected oxygen, due to leakage of air during sampling. The SCION Instruments RGA is also able to analyse natural gas. For natural gas analysis the physical properties can also be calculated in Eclipse, with extended physical property options.




Peak ID	Compound Name
1	Helium
2	Hydrogen

Figure 2a: Chromatogram of Refinery Gas Standard Front TCD Channel of RGA



Peak ID	Compound Name
1	Carbon dioxide
2	Ethylene
3	Ethane
4	Acetylene
5	Hydrogen sulfide
6	Oxygen/Argon
7	Nitrogen
8	Methane
9	Carbon monoxide

Figure 2b: Chromatogram of Refinery Gas Standard Middle TCD Channel



Peak ID	Compound Name
1	C6+
2	Methane
3	Ethane
4	Ethylene
5	Propane
6	Cyclopropane
7	Propylene
8	i-Butane
9	n-Butane
10	Propadiene
11	Acetylene
12	t-2-Butene
13	1-Butene
14	i-Butene
15	c-2-Butene
16	i-Pentane
17	n-Pentane
18	1,3-Butadiene
19	Propyne
20	Butyne

Figure 2c: Chromatogram of Refinery Gas Standard Rear FID Channel

Conclusion

SCION Instruments provides a robust, high resolution, complete solution for refinery gas streams; the Refinery Gas Analyser. This analyser determines the chemical composition of reformed gas and similar gaseous mixtures, including natural gas.

The hardware and columns of the Refinery Gas Analyser are chosen in such a way that the analysis of refinery gas and other product streams can be achieved with a standard configuration.

In combination with Eclipse, gas calculation software, an extended list of physical properties can be calculated which will give important information about the characteristics of the gas streams.

The combination of the Refinery Gas Analyser with Eclipse software is a complete solution for analysing refinery gas streams.

Peak Report

Signal 1 Front (TCD)

Time	Area	Conc	Norm	Name
3.4680	47.654e03	9.995	9.995	Helium
3.7880	41.392e03	4.949	4.949	Hydrogen

Signal 2 Middle (TCD)

Time	Area	Conc	Norm	Name
2.6863	31.455e03	0.998	0.998	Carbon Dioxide
7.5610	30.016e03	1.001	1.001	Hydrogen Sulfide
8.9417	33.100e03	0.999	0.999	Oxygen
9.8610	158.240e03	5.000	5.000	Nitrogen
12.8790	10.613e03	0.496	0.496	Carbon Monoxide

Signal 3 Rear (FID)

Time	Area	Conc	Norm	Name
3.4053	5.417e03	1.216	1.216	Further compounds with Carbon atoms ≥ 5
3.7763	23.653e03	24.960	24.960	Methane
4.2273	11.309e03	5.999	5.999	Ethane
5.0033	47.065e03	24.980	24.980	Ethylene
6.7797	14.249e03	5.017	5.017	Propane
12.9240	1.390e03	0.501	0.501	Cyclo Propane
13.2250	13.560e03	5.002	5.002	Propylene
14.9447	1.756e03	0.504	0.504	i-Butane
15.7617	3.579e03	1.004	1.004	n-Butane
16.8263	806.715e00	0.500	0.500	Propadiene
17.4103	1.717e03	1.003	1.003	Acetylene
19.5253	1.706e03	0.510	0.510	trans-2-Butene
19.8797	1.755e03	0.519	0.519	1-Butene
20.3703	3.450e03	1.005	1.005	i-Butylene
20.5983	866.598e00	0.200	0.200	neo-Pentane
20.8023	1.699e03	0.512	0.512	cis-2-Butene
21.6127	1.735e03	0.405	0.405	i-Pentane
22.1820	883.307e00	0.202	0.202	n-Pentane
23.2420	3.422e03	1.004	1.004	1,3-Butadiene
23.3843	2.150e03	1.021	1.021	Propin
27.2467	1.540e03	0.500	0.500	Butyne

100.00 100.00

Physical Properties

Metering Temperature (t2)	15	°C
Metering Pressure (p2)	101.3250	kPa
Combustion Temperature (t1)	15	°C
Combustion Pressure (p1)	101.3250	kPa
Molar Mass (dry)	26.102	g/mol
Z Mix (t2,p2, dry)	0.9951	
Carbon Content	73.67	g/100g
Rel. Ideal Density (dry)	0.90122	
Rel. Density (t,p, dry)	0.90567	
Ideal Density (dry)	1.10390	kg/m3
Real Density (t,p, dry)	1.10935	kg/m3
Inferior Cal. Value (t1, dry)	1,101.97	kJ/mol
Inferior Cal. Value (t1, dry)	4,221.83	kJ/100 g
Inferior Cal. Value (t1,V(t2,p2), dry)	46.60	MJ/m3
Real Inferior Cal. Value (t1,V(t2,p2), dry)	46.83	MJ/m3
Vapor Pressure	1,327.12	psi at 100 °F
Motor Octane Number	44.3	
Density LPG	719.5	kg/m3 at 15 °C
Abs Vapor Pressure Grade A	6,319	kPa at -10 °C
Gauge Vapor Pressure Grade A	6,218	kPa at -10 °C
Abs Vapor Pressure Grade B	6,791	kPa at -5 °C
Gauge Vapor Pressure Grade B	6,690	kPa at -5 °C
Abs Vapor Pressure Grade C	7,289	kPa at 0 °C
Gauge Vapor Pressure Grade C	7,188	kPa at 0 °C
Abs Vapor Pressure Grade D	8,363	kPa at 10 °C
Gauge Vapor Pressure Grade D	8,261	kPa at 10 °C
Abs Vapor Pressure Grade E	9,546	kPa at 20 °C
Gauge Vapor Pressure Grade E	9,445	kPa at 20 °C
Abs Vapor Pressure Grade F	12,260	kPa at 40 °C
Gauge Vapor Pressure Grade F	12,159	kPa at 40 °C
Abs Vapor Pressure 37.8	2,514	kPa at 37,8 °C
Gauge Vapor Pressure 37.8	2,413	kPa at 37,8 °C
Abs Vapor Pressure 50	3,082	kPa at 50 °C
Gauge Vapor Pressure 50	2,981	kPa at 50 °C
Carbon Emission Factor	70.46	g CO2/mol
Carbon Emission Factor	2.70	g CO2/g
Carbon Emission Factor	2.99	g CO2/m3
Carbon Emission Factor (inferior)	63.94	g CO2/MJ
Carbon Emission Factor (inferior)	63.94	ton CO2/TJ

Figure 3: Customised reporting of RGA results

Author Contact Details

Esther Van Bloois, Applications Specialist - SCION Instruments • Amundsenweg 22-24, Goes, The Netherlands



Read, Print, Share or Comment on this Article at: petro-online.com/Article

