



RAPID AND ACCURATE ANALYSIS OF SULPHUR AND CHLORINE IN BIOFUELS BY XRF

BACKGROUND

Biofuels are fuels produced from feedstocks such as vegetable oil, animal fat, used cooking oil, biomass, or a blend of these. They are not produced from fossil fuels, but rather from contemporary, human-induced processes like transesterification and hydrogenation. Biofuels are often used for blending with traditional fuels such as gasoil and gasoline. More recently, biofuels are being blended into aviation and marine fuels. However, it is possible to use certain biofuels independently without blending them with traditional fuels. It is important to note that all biofuels - blended or not - must meet certain sulphur regulatory limits. Although there are no regulations or methods for it currently, some biofuels may contain fairly high levels of chlorine which can cause corrosion damage during and after the production stages.

Biodiesels are first-generation biofuels that are made through the transesterification of vegetable oils, animal fats or used cooking oils.

Biofuels made through hydrogenation with non-food feedstocks are called second-generation biofuels or advanced biofuels.

CHALLENGE

Having the correct sulphur value for biofuels is critical as they are typically blended with fuels that have a maximum specification of 10 mg/kg for sulphur (15 mg/kg in the US). If a ULSD (Ultra Low Sulphur Diesel) is blended with a biodiesel, the maximum specification remains 10 mg/kg. Additionally, biofuel samples contain a greater concentration of oxygen than traditional fuel samples. This is important to note because oxygen absorbs XRF (X-ray Fluorescence) signals and as a result can cause analyzers to report falsely low sulphur and chlorine concentrations. For this reason, we developed a study to test real-world biofuel samples for sulphur and chlorine while using correction factors to correct for the bias caused by oxygen. To correct for oxygen content, it is important to know the actual concentration of oxygen in the sample because the correction factor is directly correlated to the amount of oxygen present in the sample.

SOLUTION

XRF delivers rapid and accurate results for testing sulphur and chlorine in biofuels, backed by international standard test methods (ASTM, ISO, etc.). Advantages of XRF technology include its non-destructive nature, easy sample preparation process, and quick results, in addition to accuracy that is on par with alternative technologies like UVF (ultraviolet fluorescence). Some relevant specifications and methods include EN 14214 (Fatty Acid Methyl Esters (FAME) for use in diesel engines and heating applications), ASTM D6751 (Standard Specification for Biodiesel Fuel Blend Stock (B100) for Middle Distillate Fuels), and ASTM D7467 (Standard Specification for Diesel Fuel Oil, Biodiesel Blend (B6 to B20)).

In this study, we will use Sindie +Cl to test sulphur and chlorine in real-world samples, including both first and second-generation biofuels as well as the traditional biofuels they are typically blended with. We will then apply correction factors to mitigate biased results. Sindie +Cl was used in this study due to its convenient ability to measure total sulphur and chlorine concurrently. The correction factors used were derived from the ASTM D7039 method, though these correction factors can still be applied when using the ISO 20884 sulphur method. As per Section 1 of ISO 20884, any sample with more than 3.7% oxygen content must be corrected for.

EXPERIMENT

For our biodiesel samples, the oxygen values have been calculated via the ester compositions of each sample. HVO does not contain oxygen as the hydrogenation process removes all functional groups, transforming the feedstock into paraffinic chains.

For this experiment the following samples were collected and analyzed for sulphur and chlorine:

- Hydrotreated Vegetable Oil (HVO)
- Rapeseed Methyl Ester (RME)
- Soybean Methyl Ester (SME)
- Tallow Methyl Ester (TME)
- Used Cooking Oil Methyl Ester (UCOME)

Each sample was measured ten times under repeatability conditions; *condition of measurement, out of a set of conditions that includes the same measurement procedure, same operator, same measuring system, same operating conditions, and same location, with replicate measurements on the same or similar objects over a short period of time.* Each sample was separated into ten aliquots via pipette into ten standard XRF cups. The samples were then sealed with Etnom® sample film, placed into Sindie +Cl, and measured for 300 seconds on a mineral oil calibration curve.

NOTE: It is possible to optimize the accuracy and precision of the results by using a custom calibration range suited to your product streams.

RESULTS

For each biodiesel sample, the oxygen content has been calculated from the ester composition, which was obtained via GC (Gas Chromatography) analysis. The results can be found in **Table 1** together with the correction factor that should be applied per element/sample. **Tables 2 and 3** display the oxygen correction factors for samples containing 0-19 wt% oxygen measured on a mineral oil calibration. Note that these factors are specific to the analyzer geometry of the Sindie® and Clora® analyzer series.

Table 1: Oxygen Content of Analyzed Samples

Sample	Oxygen Content (mass %)	Correction Factor S	Correction Factor Cl
UCOME	11.21	1.1914	1.1961
TME	11.35	1.1914	1.1961
RME	10.82	1.1914	1.1961
SME	10.91	1.1914	1.1961

Table 2: Oxygen Correction Table for Sulphur in Biofuels

Oxygen, wt%	0%	1%	2%	3%	4%	5%	6%	7%	8%	9%
0%	1.000	1.0174	1.0348	1.0522	1.0696	1.0870	1.1044	1.1218	1.1392	1.1566
10%	1.1740	1.1914	1.2088	1.2262	1.2436	1.2610	1.2784	1.2958	1.3132	1.3306

Table 3: Oxygen Correction Table for Chlorine in Biofuels

Oxygen, wt%	0%	1%	2%	3%	4%	5%	6%	7%	8%	9%
0%	1.000	1.0178	1.0357	1.0535	1.0713	1.0891	1.1070	1.1248	1.1426	1.1604
10%	1.1783	1.1961	1.2139	1.2318	1.2496	1.2674	1.2852	1.3031	1.3209	1.3387

NOTE: To determine the appropriate correction factors in Tables 1 & 2, select the row that matches the most significant figure of the oxygen concentration, and then find the column that matches the least significant figure. The intersection of the row and the column is the correction factor. The correction factor is applied by multiplying the measured result by the correction factor to obtain the oxygenated correction chlorine value. Once the samples were analyzed, data was compiled for each sample as shown in Tables 4 through 12.

Table 4: HVO Concentrations (mg/kg)

Test	Sulphur (300s)	Chlorine (300s)
1	0.31	<0.30mg/kg
2	0.47	<0.30mg/kg
3	0.33	<0.30mg/kg
4	0.40	<0.30mg/kg
5	0.49	<0.30mg/kg
6	0.40	<0.30mg/kg
7	0.47	<0.30mg/kg
8	0.29	<0.30mg/kg
9	0.51	<0.30mg/kg
10	0.43	<0.30mg/kg
average	0.41	-
st. dev.	0.08	-
RSD%	19.07	-

NOTE: The chlorine content for this particular HVO sample is very low. This sample has also been measured on Clora 2XP for 300 seconds, which resulted in a concentration of < 0.07 mg/kg

Table 5: RME Concentrations (mg/kg)

Test	Sulphur (300s)	Sulphur Corrected	Chlorine (300s)	Chlorine Corrected
1	3.34	3.98	10.25	12.26
2	3.03	3.61	13.28	15.88
3	3.35	3.99	12.35	14.77
4	3.41	4.06	10.69	12.79
5	3.34	3.98	10.30	12.32
6	3.48	4.15	11.23	13.43
7	3.44	4.10	12.37	14.80
8	3.28	3.91	12.02	14.38
9	2.94	3.50	8.60	10.29
10	3.28	3.91	9.97	11.93
average	3.29	3.92	11.11	13.28
st. dev.	0.17	0.21	1.41	1.68
RSD%	5.28	5.28	12.68	12.68

NOTE: ASTM D2622 has a sulphur scope of 3 mg/kg to 4.6 wt% of total sulphur. ISO 20884 has a sulphur scope of 5 to 500 mg/kg. The sulphur concentration of this sample is below the scope; yet Sindie +CI is still able to get reliable data. If a better precision is required at this level, our Sindie Gen 3 model has an LOD of 0.15 mg/kg in mineral oil for sulphur at 300 seconds.

Table 6: SME Concentrations (mg/kg)

Test	Sulphur (300s)	Sulphur Corrected	Chlorine (300s)	Chlorine Corrected
1	0.34	0.41	0.36	0.43
2	0.43	0.51	0.28	0.33
3	0.36	0.43	0.38	0.45
4	0.34	0.41	0.22	0.26
5	0.67	0.80	0.35	0.42
6	0.36	0.43	0.23	0.28
7	0.52	0.62	0.26	0.31
8	0.40	0.48	0.42	0.50
9	0.51	0.61	0.28	0.33
10	0.47	0.56	0.52	0.62
average	0.44	0.52	0.33	0.39
st. dev.	0.11	0.13	0.09	0.11
RSD%	23.91	23.91	28.57	28.57

Table 7: TME Concentrations (mg/kg)

Test	Sulphur (300s)	Sulphur Corrected	Chlorine (300s)	Chlorine Corrected
1	8.36	9.96	0.45	0.54
2	8.31	9.90	0.49	0.59
3	8.57	10.21	0.37	0.44
4	8.06	9.60	0.53	0.63
5	7.88	9.39	0.38	0.45
6	7.91	9.42	0.52	0.62
7	8.76	10.44	0.62	0.74
8	7.66	9.13	0.48	0.57
9	8.78	10.46	0.47	0.56
10	7.50	8.94	0.49	0.59
average	8.18	9.74	0.48	0.57
st. dev.	0.45	0.53	0.07	0.09
RSD%	5.48	5.48	15.06	15.06

Table 8: UCOME Concentrations (mg/kg)

Test	Sulphur (300s)	Sulphur Corrected	Chlorine (300s)	Chlorine Corrected
1	7.61	9.07	421	504
2	8.38	9.98	424	507
3	7.61	9.07	422	505
4	7.86	9.36	425	509
5	8.18	9.75	426	510
6	7.23	8.61	429	514
7	7.07	8.42	429	514
8	7.93	9.45	427	511
9	7.88	9.39	423	506
10	7.61	9.07	426	510
average	7.74	9.22	425	509
st. dev.	0.40	0.48	2.88	3.45
RSD%	5.16	5.16	0.68	0.68

CONCLUSION

With the continued development of biofuels and a push for higher concentrations in traditional fuels, petroleum professionals are looking to utilize technology that delivers rapid and accurate results for on-site biofuel measurement. Sindie +CI is a viable solution, delivering total sulphur and chlorine in one measurement without the need for a matrix-matched calibration by simply applying a correction factor to the results. This allows professionals to certify their biofuel products more efficiently than with other methods. Sindie +CI complies with ASTM D2622, and therefore meets the regulatory limit of ASTM D6751 and ASTM D7467. It also meets the precision requirements for ISO 20884 and can be used for EN 14214.

View the full study by visiting xos.com/biofuelsXRF

Contact Details - 1.518.880.1500 • info@xos.com • xos.com



SINDIE +CI

Sindie +CI delivers exceptional reproducibility for both sulphur and chlorine analysis with one push of a button and zero hassle. Samples are measured directly, which means it can analyze even the heaviest of hydrocarbons like crude oil or coker residuals, without the hassle of boats, injectors, furnaces, or changing detectors. Sindie +CI is powered by Monochromatic Wavelength Dispersive X-Ray Fluorescence (MWDXRF®): an elemental analysis technique offering significantly enhanced detection performance over traditional XRF technology. Using the industry's most advanced optics, doubly curved crystals, Sindie +CI achieves a high signal-to-background ratio and delivers very precise measurements of low sulphur and chlorine.

SINDIE ONLINE

Sindie® Online is an industrial grade process sulphur analyzer with breakthrough detection capability to monitor ultra-low sulphur in biofuels. This process analyzer can take a measurement every 30 seconds, enabling refiners to finetune their process control and maximize profits. Powered by MWDXRF®, Sindie Online uses ASTM D7039 technology, the same technology used in Sindie benchtop. Sindie Online offers resilience to feedstock changes and the ability to measure a variety of challenging sample types without the need for significant changes to sample conditioning.



SINDIE 2622

Sindie 2622 complies with ASTM D2622, D7039 and ISO 20884 methods, enabling complete flexibility in sulphur analysis. No compromises in detection, performance and reliability — it is the ideal sulphur analytic solution from crude feedstocks and intermediates to low sulphur biofuels.

Have questions around applications, method compliance, or best practices? Whatever you need, we're here to help. Sign up today at xos.com/virtualchats.

CHAT WITH AN XOS EXPERT, LIVE!