

# High-precision oxygen concentration analysis using a novel high-temperature pyrolysis technique

November, 17<sup>th</sup> 2015

PEFTEC 2015

Dr. Marc Ruppenthal

Product Manager Elemental Analysis  
and EA-IRMS interfaces



# Content

## 1 About ElementarGroup

2 Why oxygen?


3 rapid OXY cube features- Get OXYted!

4 Applications

5 Summary

# 110 years tradition... Heraeus

## Verbrennungsofen



System Heraeus

**Zweck des Ofens:** Ausführung von Verbrennungen bei der organischen Elementaranalyse.



# The elementargroup – tradition means trust

- Approx. 1860: W. C. Heraeus utilizes high temperature technologies in order to process platinum
- Approx. 1900: Heraeus produces high purity quartz glass which is soon used in elemental analysis




Wilhelm Carl Heraeus




Pharmacy "White Unicorn"

# The elementargroup – tradition means trust



- Approx. 1905:


First generation of elemental analyzers with electrical furnaces produced by Heraeus



- 1923:

Fritz Pregl receives the Nobel Prize in Chemistry for his invention of the *method of micro-analysis of organic substances*, **using dedicated analytical equipment manufactured by Heraeus !**

# The elementargroup history



# The elementargroup history



EXCELLENCE IN ISOTOPEs

Isoprime Ltd was acquired



vario TOC cube



trace SN cube



rapid N exceed, rapid OXY cube

2006

2009

2012

2008

2010

2014

## Introduction of the "cube" platform

vario MICRO cube, vario EL cube,  
vario MACRO cube, rapid N cube



rapid CS cube, vario ISOTOPE cube,  
vario PYRO cube



vario MAX cube, iso TOC cube



# The elementargroup history



rapid MAX N exceed, rapid MiCRO N cube

2015

## Introduction of the inductor series

inductar CS cube, inductor ONH cube,  
inductar EL cube






# The elementargroup – about us

elementargroup



- The technology leader in
  - Elemental Analysis (EA)
  - Isotope Ratio Mass Spectrometry (IRMS)
- Offers an unmatched combination of
  - innovative technologies,
  - tailor-made solutions and
  - comprehensive support.

# The elementargroup – key markets



Chemicals



Agriculture



Energy



Forensics



Materials



Environmental



# Some of our customers



# Content

1 About ElementarGroup

**2 Why oxygen?**

3 rapid OXY cube features- Get OXYted!

4 Applications

5 Summary

# Why Oxygen Analysis?



- Oxygen is an important parameter to describe:
  - Coke and coal (ASTM D3176)
  - Gasoline and methanol fuels (ASTM D5622)
  - Polymers (fire/decomposition resistance)
  - Fine chemicals (quality control)
- Two approaches:
  - Indirect:  $100\% - \text{CHNS}\% - \text{ash}\% = \text{O}\%$
  - Direct measurement

# Indirect Oxygen Analysis

- Indirect approach:  $100\% - \text{CHN}\% - \text{S}\% - \text{ash}\% = \text{O}\%$ 
  - Drawbacks:
    - Time and resources consuming
    - Summation of errors
    - Up to three instruments required

*"The result so obtained is affected by errors incurred in the other determinations of the ultimate analysis and also by the changes in weight of the ash-forming constituents of ignition."* (ASTM D3176)



# Direct Oxygen Analysis

- New applications require more precise, direct determination of oxygen

*"The presence of oxygen-containing compounds in gasoline can promote more complete combustion, which reduces carbon monoxide emissions. The Clean Air Act (1992) requires that gasoline sold within certain, specified geographical areas contain a minimum percent of oxygen by mass (presently 2.7 mass %) during certain portions of the year. The requirement can be met by blending compounds such as methyl tertiary butyl ether, ethyl tertiary butyl ether, and ethanol into the gasoline."* (ASTM D5622-95)




- Increasing need for a direct and precise determination of oxygen concentrations, especially in the solid and liquid fuel industry

# Direct Oxygen Analysis: Reductive Pyrolysis

- Sample is pyrolyzed in helium atmosphere in presence of carbon black
  - Sample oxygen (O) is converted to carbon monoxide (CO)
  - CO is detected
  - Boudouard equilibrium:  
$$\text{CO}_2 + \text{C} \rightarrow 2 \text{CO}$$
$$2 \text{CO} \rightarrow \text{CO}_2 + \text{C}$$
- Only quantitative conversion to CO results in matrix-independence



# Boudouard equilibrium



# We want CO! Simplified Pyrolysis Chemistry

- The higher the temperature, the more CO!
- Temperatures around 1200° C for sub-percent CO<sub>2</sub> formation
- But: Temperatures > 1200° C require unique technical solutions!

**Depending on temperature, a mixture of CO and CO<sub>2</sub> is produced  
→ biased oxygen measurement (matrix effects!)**

# Get OXYted about oxygen analysis - rapid OXY cube -



- Temperature 1450° C
- **Matrix independent**, no CO<sub>2</sub> formation
- No matrix standards needed
- Excellent, highly accurate results
- Ready to go, stand-alone analyzer

# Content

1 About ElementarGroup


2 Why oxygen?

**3 rapid OXY cube features- Get OXYted!**

4 Applications

5 Summary


# Get OXYted! - Blank free pyrolysis reactor -




# Get OXYted! - Patented backflush -

- **Problem:** Pyrolysis is not as well defined as combustion. Interferences can lead to wrong results.
- Function of backflush system:
  - Adsorption of CO during pyrolysis phase
  - Desorption with clean helium and temperature ramp
  - Tailing from pyrolysis byproducts does not enter the TCD but goes to vent
  - Focused peak due to temperature ramp
- **Result:** Full baseline separation of CO from interfering gases and full recovery

patented backflush




# Analysis of caffeine without backflush



No.: 89 IAEA 600



# Analysis of caffeine with backflush



**Highest accuracy  
is guaranteed!**

# Content

- 1 About ElementarGroup
- 2 Why oxygen?
- 3 rapid OXY cube features- Get OXYted!
- 4 Applications**
- 5 Summary

# Oxygen analysis has never been more reliable


SAMPLE	O [%]	ABS. SD [%]
OIL 1	0,34	0,03
OIL 2	0,66	0,04
COKE 1	3,93	0,05
COKE 2	1,88	0,10
COAL	6,32	0,08
CAFFEINE	16,94	0,07
ACETANILIDE	11,97	0,05
SULFANILIC ACID	27,77	0,07
SIVER PHOSPHATE	15,35	0,15
BARIUM SULFATE	27,72	0,44
FLUOROBENZOIC ACID	23,50	0,12

# rapid OXY cube liquid injection

- rapid OXY cube + VLS




# Liquid samples: Proof of principle



- 34 consecutive hexanol injections over wide range of injection volumes (3-9  $\mu\text{l}$ )
- Theoretical value: 15.66%
- Mean: 15.65%
- Abs. SD: 0.03%
- Excellent results
- Excellent repeatability
- Excellent accuracy

# Oxygen Content in Gasoline Blending



- Oxygen content corresponds to unwanted CO emission in engines
- Minimum oxygen content is legal requirement (e.g. Clean Air Act)
- Pure fossil gasoline with low oxygen content is blended with high oxygen content biofuels
- **Excellent repeatability**
- **No memory effects**

# Content

- 1 About ElementarGroup
- 2 Why oxygen?
- 3 rapid OXY cube features- Get OXYted!
- 4 Applications
- 5 Summary**

rapid OXY  cube

*Get OXYted about oxygen analysis*



High sensitivity




High data quality



High sample throughput



Extreme durability

rapid OXY  cube

## Highlights

- Highest accuracy and data reliability
- Patented backflush for blank-free oxygen analyses
- Designed for unattended 24/7 operation
- Unmatched operating comfort through sophisticated self-diagnosis
- High robustness – low maintenance





Good to know,  
that we are a strong family.



As a family enterprise we stand for  
highest reliability and personal commitment.



**Elementar thanks  
for your attention !**