



... efficient, environmentally friendly and cost effective approach

Single Tube Sampling: Analysis of Volatile and Semi-Volatile Organics in Air.

The Cost Effective

Green Solution



*Timon Huybrighs, PerkinElmer
Lee Marotta, PerkinElmer
Roberta Provost, Pace Analytical Services*

- ▶ Introduction
- ▶ Industry Problem/Solution
 - Combining a two analysis solution into one analysis
- ▶ EPA Method TO-17 and Thermal Desorption
- ▶ Analytical Data and Performance
- ▶ Results from Site Sampling (MGP site)
- ▶ Air Monitoring Study Summary



New Sorbent Tubes designed by PerkinElmer



- ▶ 2010: Soil Vapor Intrusion (SVI) Tube (patented)
 - (C_3 to C_{26})
 - Combines VOC & SVOC from the seven VOA gases to pyrene
 - Thank you to CARO Analytical Services for their help

- ▶ 2011: XRO-640 (patent pending)
 - (C_6 to C_{44})
 - Residue in Liquefied Petroleum Gas (LPG)
 - Combines VOC & SVOC from BTEX to benzo(g,h,i)perylene
 - Thank you to Alberta Innovates Tech Futures for their help

- ▶ 2013: XRO-440 (patent pending)
 - (C_4 to C_{44})
 - Combines VOC & SVOC from 1,3-butadiene to benzo(g,h,i)perylene
 - Thank you to Pace Analytical Services for their help



- Electronic control of all flows
 - Programmed flow, velocity or pressure
 - Enables consistent RT precision
- **Automates** spiking internal standard as a gas onto each tube
- **Automates** spiking a surrogate prior to sending tubes out for sampling
- **Automates** sample tube and cold trap impedance check to validate trap and tube
- **Automates** sample recollection - confirmatory analysis through sample recollection on the same or new tube
- **Automates** tube conditioning during analysis
- **Automates** leak check of tube and trap prior to each analysis
- Excellent water management



...ease of use, accurate, precise

PM10 (particulates)



TO-15 (VOCs)



TO-13 (PAHs)



Goal: one analysis instead of two


- Eliminate liquid extractions using TO-17
 - Save time and \$\$\$
 - Improve productivity and efficiency
 - Enhance recoveries
 - A **Greener** analysis

TO - 13



TO - 15

TO - 17



Four Canisters to a box ... Yikes!



... can fit over 500 Tubes!

The Analytical Solution for Air Monitoring



... efficient, environmentally friendly and cost effective approach

EPA Method TO-17


***to include SVOC to
benzo(g,h,i)perylene
(16 regulated PAH)***



Research required *NEW* TD Tube Design

Volatiles	Semi-volatiles
1,3-Butadiene	Napthalene
Benzene	2-Methylnapthalene
Toluene	1-Methylnapthalene
Ethyl Benzene	Acenaphthylene
Xylenes	Acenaphthene
	Fluorene
	Phenathrene
	Acenaphthene
	Fluoranthene
	Pyrene
	Benzo(a)anthracene
	Chrysene
	Benzo(b)fluoranthene
	Benzo(k)fluoranthene
	Benzo(e)pyrene
	Benzo(a)pyrene
	Indeno(1,2,3-c,d)pyrene
	Dibenz(a,h)anthracene
	Benzo(g,h,i)perylene

Example Tube



... Started with these targets for MGP sites

The Analytical Solution for Air Monitoring




... efficient, environmentally friendly and cost effective approach




Introduction to Functioning

How it works?

Step 1: Sample Tube (Primary) Desorption



Step 2: Transfer Sample to Analytical Column







... efficient, environmentally friendly and cost effective approach

Note: Data was acquired in Simultaneous Full Scan/ SIM Mode. Only results from Full Scan are reported

Results of Analytical Performance



Total Ion Chromatogram (TIC)




Calibration, Precision and Reporting Limit

Target Compound	Range 0.2 to 50ng	Reporting Limit (ug/m ³) 45L sample volume	Precision (%RSD) n=6
1,3-Butadiene	0.9961	0.0111*	1.89
Benzene	0.9971	0.0044	0.90
Toluene	0.9991	0.0044	0.94
Ethyl Benzene	0.9989	0.0044	0.77
m & p - Xylenes	15.54%	0.0044	0.95
o - Xylene	0.9994	0.0044	1.57
Naphthalene	25.07%	0.0044	0.92
2-Methylnaphthalene	11.79%	0.0044	1.69
1-Methylnaphthalene	19.05%	0.0044	0.65
Acenaphthylene	11.32%	0.0044	1.87
Acenaphthene	14.40%	0.0044	1.48
Fluorene	20.96%	0.0044	2.27
Phenanthrene	8.13%	0.0044	1.67
Anthracene	15.54%	0.0044	2.27
Fluroanthene	7.23%	0.0044	1.41
Pyrene	22.44%	0.0044	1.24
Benzo[a]anthracene	18.93%	0.0044	2.04
Chrysene	19.21%	0.0044	1.92
Benzo[b&k]fluoranthene	16.21%	0.0044	5.96
Benzo[e]pyrene	16.61%	0.0044	0.80
Benzo[a]pyrene	10.86%	0.0044	0.99
Indeno[1,2,3-c,d]pyrene	20.28%	0.0044	1.78
Dibenz[a,h]anthracene	0.9951	0.0044	1.21
Benzo[g,h,i]perylene	0.9952	0.0044	1.97



o-Xylene


Compound 12 name: o - Xylene
Coefficient of Determination: 0.999930
Calibration curve: $0.292679 * x + 0.00222540$
Response type: Internal Std (Ref 1), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



Pyrene

Compound 22 name: Pyrene
Coefficient of Determination: 0.999243
Calibration curve: $1.90286 * x + -0.239747$
Response type: Internal Std (Ref 4), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None

Correlation Coefficient
 $R^2 = 0.9992$



The Analytical Solution for Air Monitoring



... efficient, environmentally friendly and cost effective approach




***Breakthrough and Recovery
Experiments and Results***

What is Breakthrough?

- ▶ Occurs when target compounds are not adsorbed by adsorbents
- ▶ EPA TO-17 definition: "The volume sampled when the amount of analyte collected in a back-up sorbent tube reaches a certain percentage (typically 5%) of the total amount collected by both sorbent tubes"



First Breakthrough Experiments



- 8000ng on tube
- 85% humidified N₂ for 100 min at a flow of 100 mL/min
- Total volume 10 L sampling
- All studies performed in triplicate



2nd Breakthrough Experiment and Results

- ▶ A primary TD tube was attached to a gaseous standard to continuously deliver target compounds (mimics a real-world sampling event)
- ▶ A BT tube was attached and monitored on a regular basis
- ▶ Ultimately, the primary tube was loaded with >200mg analyte with no detectable breakthrough

Target Analyte	% BT	% BT	% BT
1,3-Butadiene	nd	nd	nd
Benzene	nd	nd	nd
Toluene	nd	nd	nd
Ethyl Benzene	nd	nd	nd
m & p - Xylenes	nd	nd	nd
o - Xylene	nd	nd	nd
Naphthalene	nd	nd	nd
2-Methylnaphthalene	nd	nd	nd
1-Methylnaphthalene	nd	nd	nd
Acenaphthylene	nd	nd	nd
Acenaphthene	nd	nd	nd
Fluorene	nd	nd	nd
Phenanthrene	nd	nd	nd
Anthracene	nd	nd	nd
Fluoroanthene	nd	nd	nd
Pyrene	nd	nd	nd
Benzo[a]anthracene	nd	nd	nd
Chrysene	nd	nd	nd
Benzo[b&k]fluoranthene	nd	nd	nd
Benzo[e]pyrene	nd	nd	nd
Benzo[a]pyrene	nd	nd	nd
Indeno[1,2,3-c,d]pyrene	nd	nd	nd
Dibenz[a,h]anthracene	nd	nd	nd
Benzo[g,h,i]perylene	nd	nd	nd




Carryover and Recovery

Recovery/Carryover Experiments

- Analyzed spiked tube (50ng)
- Analyzed trap
- Analyzed valve
- Re-analyzed spiked tube

Target Analyte	Trap Test	Tube Test	Valve Test
Benzene	nd	nd	nd
Toluene	nd	nd	nd
Ethylbenzene	nd	nd	nd
m&p-Xylene	nd	nd	nd
o-Xylene	nd	nd	nd
Naphthalene	nd	nd	nd
2-Methylnaphthalene	nd	nd	nd
Acenaphthylene	nd	nd	nd
Acenaphthene	nd	nd	nd
Fluorene	nd	nd	nd
Phenanthrene	nd	nd	nd
Fluoranthene	nd	nd	nd
Chrysene	nd	nd	nd
Benzo[a]pyrene	nd	nd	nd
Indeno[1,2,3-cd]pyrene	nd	nd	nd
Benzo[g,h,i]perylene	nd	nd	nd



The Analytical Solution for Air Monitoring



... efficient, environmentally friendly and cost effective approach




Site Study

TO-17 data was collected in simultaneous Full Scan/SIM mode
only Full Scan data is presented

- ▶ Compare TO-13 / TO-15 to TO-17: results from an active MGP remediation site
- ▶ 72-hour sample collection
- ▶ Continuous sampling for six weeks
- ▶ Two sample locations selected (AMS-01 and AMS-03)
- ▶ Three 72-hour samples from each site were selected for comparison



Site Map



Site Setup



TO-17 Site Study expanded view

- ▶ Two types of tubes investigated (XRO-444 and XRO-644)
- ▶ Each type was sampled in duplicate
- ▶ One of the duplicates had filter attached and analyzed
- ▶ A breakthrough tube was attached to every tube sampled



72 Hour Sampling Duration

- ▶ TO-13 = ~1,000,000 Liters (1000 m³)
- ▶ TO-15 = 6 Liters
- ▶ TO-17 = ~45 Liters



Moisture: Hydrophobic adsorbents

Tube	Time for Dry Purge
Sample Tube Type 1	nd added moisture on tube
Sample Tube Type 2	2 minute Dry Purge

Only slight water retention with
45L sample volume!




- ▶ There was non-detectable targets on the breakthrough tubes from the site studies with an average of 45 liter sample volume



Duplicate Concentrations and On-Tube Values

Target Analyte	ug/m3 (first tube)	ug/m3 (second tube)	equates to ng on tube	% relative dif
1,3-Butadiene	0.0193	0.0200	0.87	3.6
Benzene	0.3827	0.4090	17.2	6.6
Toluene	1.3329	9.7246	60.0	152
Ethyl Benzene	0.1543	0.2136	6.94	32.2
m & p - Xylenes	0.4413	0.6081	19.9	31.8
o - Xylene	0.1434	0.1586	6.45	10.1
Naphthalene	3.1182	3.4084	140	8.9
2-Methylnaphthalene	0.6185	0.6083	27.8	1.7
1-Methylnaphthalene	0.2647	0.3138	11.9	17.0
Acenaphthylene	0.0656	0.0492	2.95	28.6
Acenaphthene	0.3022	0.2660	13.6	12.8
Fluorene	0.1238	0.1768	5.57	35.3
Phenanthrene	0.0547	0.0931	2.46	52.0
Anthracene	0.0803	0.0915	3.61	13.1
Fluroanthene	0.0040	0.0050	0.18	21.7
Pyrene	0.0032	0.0050	0.14	45.4
Benzo[a]anthracene	0.0067	0.0054	0.30	21.0
Chrysene	0.0046	0.0033	0.21	32.7
Benzo[b&k]fluoranthene	0.0044	nd	0.20	
Benzo[e]pyrene	0.0044	nd	0.20	
Benzo[a]pyrene	0.0074	nd	0.33	
Indeno[1,2,3-c,d]pyrene	nd	nd		
Dibenz[a,h]anthracene	nd	nd		
Benzo[g,h,i]perylene	nd	nd		




Results from the Study – AMS-01



Sample ID	AMS-01-081313		AMS-01-081613		AMS-01-082213	
	TO13 & 15	TO17	TO13 & 15	TO17	TO13 & 15	TO17
Benzene	0.68*	0.38	0.50	0.56	1.9	0.75
Ethyl Benzene	1.8	0.15	1.4	0.29	5.1	0.66
Toluene	15	1.3	1.6	1.1	8.3	1.6
m,p-Xylene	3.7	0.44	2.8	0.83	10	1.53
o-Xylene	1.80	0.14	1.4	0.24	5.1	0.47
Naphthalene (TO15)	1.4	3.1	2.9	1.1	6.2	1.10
Naphthalene (TO13)	0.68	3.1	0.68	1.1	0.082	1.10
2-Methylnaphthalene	0.25	0.62	0.33	0.26	0.044	0.26
1-Methylnaphthalene	0.12	0.26	0.17	0.16	0.031	0.29
Acenaphthylene	0.0058	0.066	0.0073	0.045	0.011	0.065
Acenaphthene	0.13	0.30	0.18	0.14	0.039	0.42
Fluorene	0.070	0.12	0.10	0.042	0.037	0.16
Phenanthrene	0.076	0.055	0.11	0.015	0.065	0.11
Anthracene	0.0039	0.080	0.0039	0.0040	0.0049	0.17
Fluoranthene	0.0092	0.0040	0.014	0.0047	0.021	0.0076
Pyrene	0.0050	0.0032	0.0073	0.0047	0.010	0.0044
Benzo(a)anthracene	0.0006	0.0067	0.0013	0.0047	0.00051	0.0044
Chrysene	0.00089	0.0046	0.0017	0.0047	0.0010	0.0044
Benzo(b+k)fluoranthene	0.00092	0.0088	0.0033	0.014	0.00062	0.0094
Benzo(e)pyrene	0.00051	0.0044	0.0013	0.0084	0.00050	0.0044
Benzo(a)pyrene	0.00048	0.0074	0.0013	0.0047	0.00050	0.0044
Indeno(1,2,3-cd)pyrene	0.00046	0.0044	0.0011	0.0047	0.00050	0.0044
Dibenz(a,h)anthracene	0.00046	0.0044	0.00050	0.0075	0.00050	0.0044
Benzo(g,h,i)perylene	0.00050	0.0081	0.0016	0.0047	0.00050	0.0044


NOTE: yellow cells are non-detect with the reporting limit value for that target in the cell

Chart for Site AMS-01-081313




Note: dotted, unfilled bars are non-detects. Value represents reporting limit.

Chart for Site AMS-01-081613



Note: dotted, unfilled bars are non-detects. Value represents reporting limit.

Chart for Site AMS-01-082213



Note: dotted, unfilled bars are non-detects. Value represents reporting limit.


Results from the Study – AMS-03



Sample ID Analyte	AMS-03-081313		AMS-03-081613		AMS-03-082213	
	TO13 & TO15	TO17	TO13 & 15	TO17	TO13 & 15	TO17
Benzene	0.52	0.54	1.20	0.56	1.2	0.48
Ethyl Benzene	1.4	0.29	1.4	0.29	3.2	0.78
Toluene	2	1.9	1.6	1.1	4.4	0.86
m,p-Xylene	2.8	0.87	2.8	0.83	6.5	0.91
o-Xylene	1.4	0.22	1.4	0.24	3.2	0.26
Naphthalene (TO15)	3.2	6.9	3.3	1.1	3.9	2.4
Naphthalene (TO13)	1.2	6.9	0.95	1.3	0.18	2.4
2-Methylnaphthalene	0.48	0.74	0.51	0.41	0.10	0.67
1-Methylnaphthalene	0.25	0.36	0.26	0.25	0.088	1.2
Acenaphthylene	0.023	0.12	0.0094	0.047	0.17	0.40
Acenaphthene	0.27	0.25	0.30	0.21	0.17	1.7
Fluorene	0.14	0.10	0.17	0.081	0.13	0.45
Phenanthrene	0.16	0.077	0.17	0.019	0.24	0.24
Anthracene	0.023	0.11	0.0061	0.0258	0.016	0.36
Fluoranthene	0.023	0.0092	0.017	0.0047	0.13	0.013
Pyrene	0.023	0.0043	0.0088	0.0047	0.027	0.0018
Benzo(a)anthracene	0.023	0.0051	0.0011	0.0047	0.00074	0.00029
Chrysene	0.023	0.0051	0.0016	0.0047	0.0014	0.0046
Benzo(b+k)fluoranthene	0.046	0.010	0.0027	0.0094	0.00091	0.0092
Benzo(e)pyrene	0.023	0.0051	0.0011	0.0047	0.00050	0.0046
Benzo(a)pyrene	0.023	0.0051	0.0010	0.0066	0.00050	0.0046
Indeno(1,2,3-cd)pyrene	0.023	0.0051	0.0009	0.0065	0.00050	0.0046
Dibenz(a,h)anthracene	0.023	0.0051	0.00050	0.0140	0.00050	0.0046
Benzo(g,h,i)perylene	0.023	0.0051	0.0013	0.0220	0.00050	0.015


NOTE: yellow cells are non-detect with the reporting limit value for that target in the cell

Chart for Site AMS-03-081313




Note: dotted, unfilled bars are non-detects. Value represents reporting limit.

Chart for Site AMS-03-081613



Note: dotted, unfilled bars are non-detects. Value represents reporting limit.

Chart for Site AMS-03-082213



Note: dotted, unfilled bars are non-detects. Value represents reporting limit.

- ▶ EPA Method TO-17 is performance-based, guidance method
 - Section 2.5 states: "...This method provides performance criteria to demonstrate acceptable performance of the method (*or modifications of the method*) for monitoring a compound or set of compounds."

- ▶ EPA has seen this data and has given verbal acceptance stating that TO-17 is performance based so targets may be included as long as criteria is met
 - U.S. Environmental Protection Agency
Office of Air Quality Planning and Standards
Ambient Air Monitoring Group C304-06
Research Triangle Park, NC 27711



- ▶ Analytical performance proves concept
- ▶ Site data suggests this is a better alternative
- ▶ One analysis instead of two:
 - Reduce sampling and analytical costs and disposal
 - Save on shipping and labor costs
 - Enhance productivity and efficiency
 - Increase profits
 - Better for our environment ... **A Greener analysis**
- ▶ More data is available



Acknowledgements



- ▶ Amy Jacobson, Specialty Analytical Services Manager, Pace Analytical Services
- ▶ Nathan Eklund, Program Manager – Specialty Analytical Services, Pace Analytical Services
- ▶ Mariah Peronto, Account Manager, Pace Analytical Services
- ▶ James Day, Service Engineer, PerkinElmer



The Analytical Solution for Air Monitoring



... efficient, environmentally friendly and cost effective approach



Thank You !

Questions please ????

Timon.Huybrighs@perkinelmer.com

Lee.Marotta@perkinelmer.com